

# General Equipment



# Front Distribution

Medical Equipment

[info@frontdistribution.com](mailto:info@frontdistribution.com)

[www.frontdistribution.com](http://www.frontdistribution.com)

## Manual blood pressure monitor

CODE	Size	Pack Qty
101119	12 x 35 cm	1
101120	15 x 43 cm	1



## Manual Heine blood pressure monitor G5

- Micro-filter that protects the valve and manometer
- Large and clear scale, 56mm in diameter
- Calibrated
- Rapid inflation of the cuff
- Accurate screw valve for precise air emissions

CODE	Size	Pack Qty
440120	8 x 24 cm (children)	1
440132	12 x 28 cm (adult)	1



## Blood pressure monitor Omron M7 Intelli IT

- OMRON M7 Intelli IT has a mobile application that connects to Omron Bluetooth devices and provides an easy view of your measurements, history, graphs and averages.
- It is equipped with Intelli Wrap cuff.
- The cuff guarantees a much more reliable result in any position around the arm (360 degrees) which reduces the risk of incorrect results.
- The cuff is preformed and conforms comfortably to the arm creating an even pressure. It fits most arms, up to 42 cm in circumference

CODE	Pack Qty
816085	1



## Blood pressure monitor Omron MIT5s Connect Professional

- Blood pressure monitor MIT5s measures blood pressure using an oscillometric method, which means that the monitor detects blood flow through the artery and converts the movements into a digital reading. With an oscillometric monitor, a stethoscope is not needed so the monitor is easy to use.
- MIT5s combines advanced Omron technology with modern design. Its unique combination guarantees accurate measurement and allows you to track the blood pressure result on your smartphone.



CODE	Pack Qty
120223	1

## Blood pressure monitor Omron HBP-1100 (for hospitals)

- Large and clear display
- Wearable cuff in size M (22-32 cm) included
- Quality indicator of the results - if the person has moved a lot during the measurement, inflation stops and if the person has an irregular heart rate (shown with pictures).
- HBP has a zero indicator that ensures that the measurement always starts from zero
- Clinically validated in accordance with International Protocols



## Inhaler Pari Boy Classic

- All-round for acute and chronic respiratory diseases from 4 years old.
- PARI BOY Classic can be used to treat conditions such as bronchial asthma, allergic asthma and acute or chronic bronchitis.
- Inhalation times of less than 5 minutes with 2.5 ml of inhalation solution.
- PARI LC SPRINT Nebulizer for effective lung deposition even with short inhalation times.



CODE	Pack Qty
106020	1

## Maxin air Nebulizer (Complete package)

- The Maxin package contains a complete Maxin nebulizer with hose, quick coupling, cleaning needle, five nozzles and 400 dose sachets.

CODE	Pack Qty
3100	1



## Maxin mask for nebulizer

- Face-shaped, soft silicone.

CODE	Size	Pack Qty
3600	80 x 105 mm (children)	1
3601	90 x 115 mm (adults)	1



## AccuTrend Plus Kit

- Fast and easy patient-centered system for quantitative determination of three cardiovascular risk factors: glucose, cholesterol and triglycerides.
- Through a simple stick in the finger (capillary blood), the analyzes can be obtained within a few minutes.
- The instrument can also measure lactate, which makes it useful in e.g. sports medicine.
- The kit contains:
  - 1 x AccuTrend Plus device
  - 1 x Protective case for AccuTrend Kit and Test Strips
  - 1 x 200 Safe-T-Pro Plus launcher
  - 1 x 25 test strips for AccuTrend triglycerides
  - 1 x 25 AccuTrend Cholesterol Test Strips
  - 1 x 25 BM-lactate test strips
  - 1 x 50 BM AccuTrend Glucose Test Strips
  - 2 x 1.5 ml AccuTrend Cholesterol Control Solution
  - 2 x 4 ml BM lactate control solution
  - 2 x 4 ml AccuTrend glucose control solution
  - All delivered test strips have at least five months before expiry date.


CODE	Pack Qty
780040	1



## Glucose meter Contour XT

- The blood glucose meter is ready to use, directly from the package.
- The Contour XT meter is only used with the new Contour Next test strips.
- Contour XT is intended for people who test their blood sugar themselves, as well as for healthcare professionals who measure patients' blood sugar.

CODE	Pack Qty
121520	1



## Contour Next Glucose Test Stick

- Glucose stick for Contour XT and Contour next One.

CODE	Pack Qty
121521	50 strips/pack



## Otoscope / Ophthalmoscope set Heine mini 3000

CODE	Info	Pack Qty
440051	-	1
440054	With LED light	1



## KaWe Stethoscope Plano

- Hands-free stethoscope that fits well under the blood pressure cuff.
- Available in several different colors.

CODE	Color	Pack Qty
113714	Red	1
113715	Black	1
113716	Blue	1
113717	Green	1
113718	Yellow	1
113719	Purple	1
113720	Turquoise	1
113721	Pink	1



©KaWe

## KaWe Stethoscope Prestige (for children)

CODE	Color	Pack Qty
113713	Black	1



## KaWe Stethoscope Double Hose

CODE	Color	Pack Qty
113711	Black	1



©KaWe

## Littmann Stethoscope Electronic 3200

- Amplifying sound
- Friction noise reduction function
- Record multiple sounds from the same patient and compare them or save up to twelve 30-second sounds on the stethoscope from several different patients and put them in electronic files as needed.

CODE	Color	Pack Qty
100332	Navy blue	1
100330	Black	1
100331	Red	1



## Gentle Ear Thermometer Temp 520

- Infrared ear thermometer in ergonomic design that measures heat in 1 second
- Gentle Temp 520 can measure in both Celsius and Fahrenheit
- Stores the last 9 measurements
- Automatic shut-off 1 min after measurement
- 21 protective sleeves included

CODE	Pack Qty
816072	1



## Protective sleeves for Gentle Ear Thermometer Temp 520

CODE	Pack Qty
816073	40



## Thermometer Flex Temp Smart

- Flex Temp Smart is our thermometer with extra soft and flexible tip for better comfort.
- It measures rectally in just 10 seconds and 20 seconds when using orally.

CODE	Pack Qty
816076	1

