

Gowns and apparel

Masks

Gloves

Protective Clothing



Front Distribution

Medical Equipment

info@frontdistribution.com

www.frontdistribution.com

Face mask 3M Aura FFP2 med ventil

- Face mask class FFP2 (NPF 12)
- Flexible metal strip in the nose area allows for a better fit
- Benefit from 2 elasticated head straps
- Available in one size fits all

CODE	Info	Pack Qty
100101	FFP2 (NPF 12)	10



Face mask 3M Aura FFP2 utan ventil

- Face mask class FFP2 (NPF 12)
- Flexible metal strip in the nose area allows for a better fit
- Benefit from 2 elasticated head straps
- Available in one size fits all

CODE	Info	Pack Qty
100103	FFP2 (NPF 12)	20



Face mask 3M Aura FFP3 utan ventil

- Face mask class FFP3 (NPF 50)
- Flexible metal strip in the nose area allows for a better fit
- Benefit from 2 elasticated head straps
- Available in one size fits all

CODE	Info	Pack Qty
100106	FFP3 (NPF 50)	1



Face mask 3M class II blue 1810F

- Flexible metal strip in the nose area allows for a better fit
- Horizontal head straps for tying
- Available in one size fits all

CODE	Info	Pack Qty
100795	Horizontal head straps for tying	100



Face mask activaMask class IIR Blue

- Flexible metal strip in the nose area allows for a better fit
- 3 layer face mask class IIR
- Elastic head straps
- Available in one size fits all
- CE EN14683: 2019 - Type IIR.
- 93/42 EEC Klass 1 MDD.
- Outer layer (colored): polypropylene fleece - hydrophobic
- Middle layer: polypropylene meltblown fleece - BFE ≥98%
- Inner layer (white): polypropylene fleece - hydrophilic
- Nose clip: wire with polyethylene coating
- Ear straps: polyester / elastane

CODE	Info	Pack Qty
840021	3 layer face mask class IIR	50



Cotton glove 30 cm

- Bleached white material helps to promote a clinical and hygienic appearance
- Soft cotton ensures comfort for users with allergies, sensitive skin or eczema.
- 100% Cotton

CODE	Info	Pack Qty
700060	Small	20
700061	Medium	20
700062	Large	20



Latex gloves PF Semperguard Innercoated

- Ideal for working with fluids, instruments or wet glassware

CODE	Info	Pack Qty
434224	Small	100
434225	Medium	100
434226	Large	100
434227	X-Large	90



Latex gloves powder free

- Medical examination gloves

CODE	Info	Pack Qty
101560	X-Small	100
101561	Small	100
101562	Medium	100
101563	Large	100
101564	X-Large	100



Latex gloves with powder

- Ideal for working with fluids, instruments or wet glassware

CODE	Info	Pack Qty
434232	Medium	100
434233	Large	100
434234	X-Large	90



Nitrile gloves Latex free

- Medical examination gloves
- Latex free
- Extra elastic
- Textured fingertips for better grip

CODE	Info	Pack Qty
101310	X-Small	150
101311	Small	150
101312	Medium	150
101313	Large	150
101314	X-Large	150



Nitrile gloves Latex free

- Medical examination gloves
- Latex free
- High quality

CODE	Info	Pack Qty
101315	X-Small	100
101316	Small	100
101317	Medium	100
101318	Large	100
101319	X-Large	100



Clean Air Suit Pants Blue

- BARRIER® Clean air suit
- Helps create ultra clean air
- Absorbing soft inner layer for extra comfort
- Large pockets

CODE	Info	Pack Qty
550079	X-Small	28
550080	Small	26
550081	Medium	22
550082	Large	20
550083	X-Large	18



Clean Air Suit Shirt Blue

- BARRIER® Clean air suit
- Helps create ultra clean air
- Absorbing soft inner layer for extra comfort
- Large pockets.

CODE	Info	Pack Qty
550089	X-Small	28
550090	Small	26
550091	Medium	22
550092	Large	20
550093	X-Large	18



Personal/Patient clothes Pants + Shirt

- Unisex design
- Tough material

CODE	Info	Pack Qty
240020	Small	30
240021	Medium	30
240022	Large	30
240023	X-Large	30
240024	XX-Large	30

