

Injection Systems

Syringes

Needles



Front Distribution

Medical Equipment

info@frontdistribution.com


www.frontdistribution.com

Syringe 3-Part BD Plastipak Luer

- Used for general purpose injection and aspiration of fluids
- Barrel has clear and precise markings for accurate dosage
- Retaining ring prevents accidental plunger rod withdrawal
- Syringes are latex free elastomer to eliminate the risk of allergic reactions
- Comply with International Standards: CE marked and ISO 7886-1



CODE	Volume	Info	Pack Qty
1001	1 ml	Centered	120
1002	20 ml	Eccentric	60
1003	30 ml	Eccentric	60
1004	50 ml	Eccentric	60



Syringe 3-Part BD Plastipak Luer-Lock

- Ideal for applications where avoiding potential needle disengagement, medication leakage or tubing disconnect is crucial
- Integral Luer-Lock hub secures connection to needles or IV ports
- Barrel has clear and precise markings for accurate dosage
- Latex free rubber and contain medical grade silicone oil lubricant
- Compatible with our Sterile Needles
- Comply with International Standards: CE marked and ISO 7886-1



CODE	Volume	Pack Qty
130470	1 ml	100
130471	3 ml	200
130472	5 ml	125
130473	10 ml	100
130474	20 ml	60
130476	30 ml	60
130477	50 ml	60

Syringe 3-Part BD Emerald Luer

- Reduction in material to minimized waste disposal
- Designed for a perfect fit when connected with a needle hub or catheter, ensuring minimal discomfort to the patient
- Green stopper and bold scale markings for accurate dosage
- Flange size and shape provides stability and comfort during use
- Robust retaining ring prevents accidental plunger rod withdrawal
- Textured thumb press reduces slippage during administration
- Comply with International Standards: CE marked, ISO 7886-1 and ISO 594
- Does not contain natural rubber latex



CODE	Volume	Info	Pack Qty
1005	2 ml	Centered	100
1006	3 ml	Centered	100
1007	5 ml	Centered	100
100	10 ml	Centered	100

Oral syringe BD

- Special oral syringe where the end-piece prevents mounting of hypodermic needle
- Supplied with leak-tight cap to prevent loss of the medication contained in the syringe after preparation
- Ideal for use in pediatric and geriatric departments
- Syringes are latex free rubber and contain medical grade silicone oil lubricant
- Syringe varieties with dark barrel to protect sensitive fluids from UV-light



CODE	Volume	Info	Pack Qty
130310	1 ml	Amber-colored	100
130311	3 ml	Amber-colored	100
130312	5 ml	Amber-colored	100
130313	10 ml	Amber-colored	100
130320	1 ml	Transparent	100
130321	3 ml	Transparent	100
130322	5 ml	Transparent	100
130323	10 ml	Transparent	100

BD Microfine insulin syringe

- Microfine insulin syringe
- Syringes have a one-unit increment allowing you to deliver the correct dose more accurately
- Syringes are lubricated improving injecting and minimizing patient discomfort
- Needle is designed with a sharp edge that can penetrate the skin easily and requiring less force when pushing the injection plunger
- Needle also has a thin wall technology which can improve insulin flow rate



CODE	Volume	Info	Pack Qty
130300	0,3 ml	(30 G) 0,30 mm x 8 mm	100
130301	0,5 ml	(30 G) 0,30 mm x 8 mm	100
130302	0,5 ml	(29 G) 0,33 mm x 12,7 mm	100
130303	1 ml	(29 G) 0,33 mm x 12,7 mm	100

1 ml syringe with combined needle

- Low waste space reducing medication waste
- Lower risk of air bubbles formation
- No risk of the needle detaching the syringe
- Needle is retracted into the syringe barrel to prevent needlestick injuries



CODE	Volume	Info	Pack Qty
130432	1 ml	0,40 x 10mm	100
130431	1 ml	0,45 x 13mm	100

Syringe 3-Part B.Braun Omnifix Solo

- Hypodermic syringes for single use without needle
- Available in sizes from 2- 50 ml
- Highly transparent barrel
- High contrast black graduation
- Safe plunger backstop
- Piston with double sealing ring for slow aspiration and injection
- Does not contain natural rubber latex according to ISO 7886-1



CODE	Volume	Pack Qty
121190	2 ml	100
121191	5 ml	100
121192	10 ml	100
121193	20 ml	100
121194	50 ml	100

Syringe 3-Part B.Braun Omnifix Luer-Lock Solo

- Hypodermic syringes for single use with Luer-Lock needle fixation
- Available in sizes from 1- 50 ml
- Highly transparent barrel
- High contrast black graduation
- Safe plunger backstop
- Piston with double sealing ring for slow aspiration and injection
- Does not contain natural rubber latex according to ISO 7886-1

CODE	Volume	Pack Qty
121201	1 ml	100
121202	2 ml	100
121203	5 ml	100
121204	10 ml	100
121205	20 ml	100
121206	30 ml	100
121207	50 ml	100



Syringe B.Braun Omnifix Solo with centered catheter tip

- Transparent barrel
- High contrast graduation
- Permanent marking
- Easily readable
- Centric catheter fitting cone for special nutrition and irrigation
- Safe plunger backstop
- Piston with double sealing ring for slow aspiration and injection
- Does not contain natural rubber latex according to ISO 7886-1

CODE	Volume	Pack Qty
121195	50 ml	100



Syringe B.Braun Perfusor

- Transparent barrel
- High contrast graduation
- Permanent marking
- Easily readable
- Hypodermic syringes for single use with Luer-Lock needle fixation
- Safe plunger backstop
- Piston with double sealing ring for slow aspiration and injection
- Does not contain natural rubber latex according to ISO 7886-1
- Can be ordered with aspiration needle

CODE	Volume	Pack Qty
121196	20 ml	50



Syringe BD Plastipak with centered catheter tip

- Can be used for connection to non-ISO compatible Luer connections such as nasogastric tubes
- Barrel has clear and precise markings for accurate dosage
- Retaining ring prevents accidental plunger rod withdrawal
- Individually wrapped, single use, sterile syringes for patient safety
- Latex free elastomer and contain medical grade silicone oil lubricant (less than 0.25 mg/cm²)
- Tested during manufacture in compliance with ISO 7886-1
- The 100 ml syringe can used with Luer connections with aid of an included adapter



CODE	Volume	Pack Qty
130479	100 ml	25
130478	50 ml	60

Syringe Omnifix F-Luer Tuberkulin 1 ml

- Fine dosage syringes for single use without needle
- Highly transparent barrel
- High contrast black graduation
- Safe plunger backstop
- Piston with double sealing ring for slow aspiration and injection
- Does not contain natural rubber latex according to ISO 7886-1



CODE	Volume	Pack Qty
121180	1 ml	100

Syringe BD PosiFlush NaCl, Pre-filled Flushing Syringes

- Ready-to-use, sterile pre-filled flushing devices containing sodium chloride (isotonic NaCl 0,9 % = 9 mg/ml)
- Saline syringes suitable for use on a sterile fluid path, use for general IV flush procedures
- BD PosiFlush™ saline syringes are designed to eliminate syringe-induced blood reflux enhancing catheter maintenance protocols
- Special design of syringe barrel means the smaller 3ml syringes have the same diameter as the 10ml syringes, assuring a standardized pressure during the flushing process
- BD Luer-Lock™ tip for secure connection to the catheter, minimizing patient discomfort
- Do not contain: Latex, BPA, PVC or DEHP
- Sterilized container



CODE	Volume	Pack Qty
130480	3 ml	30
130481	5 ml	30
130482	10 ml	30

Acupuncture needles (without tube)

- The HEGU Xeno needle body is manufactured in exceptional quality Japanese surgical stainless steel.
- Micro polishing and ultrasound cleaning are part of the production process.
- In the final stage the needles are sterilized using ethylene oxide gas.
- 3 years sterilization time from manufacturing date



CODE	Dimension	Pack Qty
503073	0,20 x 15 mm	100
503074	0,25 x 15 mm	100
503075	0,30 x 15 mm	100
503076	0,30 x 25 mm	100
503077	0,30 x 30 mm	100
503078	0,30 x 40 mm	100
503079	0,30 x 50 mm	100

BD Pen needle

- Compatible with leading diabetes medication pens
- Thin-Wall technology design to make the needle wall thinner.
- The Pen needles have the same outer diameter but with an increase inner diameter for a better insulin flow
- The Pen needles are intended for single use



CODE	Dimension	Pack Qty
130330	4 mm 0,23 mm (32G)	100
130331	5 mm 0,25 mm (31G)	100
130334	6 mm 0,25 mm (31G)	100
130332	8 mm 0,25 mm (31G)	100
130333	12,7 mm 0,33 mm (29G)	100

BD Vacutainer Eclipse Signal Blood collection needle

- Intuitive design for ease of use, comfort and complete user control
- Built-in safety shield in alignment with the bevel of the needle
- Integrated holder and needle reduce risk of needlestick injuries
- Flashback chamber provides immediate confirmation of venous access
- Single-handed safety shield activation eliminates need for hard surface activation



CODE	Dimension	Info	Pack Qty
130490	0,7 x 25 mm	Black	50
130491	0,8 x 25 mm	Green	50

BD Vacutainer Eclipse Signal Blood collection needle with integrated holder

- Intuitive design for ease of use, comfort and complete user control
- Built-in safety shield in alignment with the bevel of the needle
- Integrated holder and needle reduce risk of needlestick injuries
- Flashback chamber provides immediate confirmation of venous access
- Single-handed safety shield activation eliminates need for hard surface activation



CODE	Dimension	Info	Pack Qty
130492	0,7 x 25 mm	Black	50
130493	0,8 x 25 mm	Green	50

Biopsy Needle BD Chiva

- Lightweight, thin-wall 22G needles are uniformly sharp for consistent, smooth entry
- The translucent hub provides sensitivity for accurate needle placement
- Supplied in peel-off blister pack
- Sterile
- Available in two different length



CODE	Dimension	Pack Qty
130010	0,7 x 153 mm	10
130011	0,7 x 203 mm	10

Infusion needle Introcan Safety 3

- Peripheral IV catheter with a passive fully automatic needlestick protection
- Multi-access blood control septum that helps to protect clinicians and patients from blood exposure
- It is a single-use device to generate intravascular and tissue access to sample blood, monitor blood pressure, or administer fluids and blood intravascularly
- It is also indicated for subcutaneous infusion therapies and 300 psi power injector applications

CODE	Dimension	Info	Pack Qty
121030	24G 0,7 x 19 mm	Yellow	50
121031	22G 0,9 x 25 mm	Blue	50
121032	20G 1,1 x 25 mm	Pink	50
121033	20G 1,1 x 32 mm	Pink	50
121034	18G 1,3 x 32 mm	Green	50
121035	18G 1,3 x 45 mm	Green	50



Infusion needle Vasofix safety IV Catheter

- Easy puncture at a wide range of angles
- Minimum effort of catheter insertion
- Self-activating Safety Shield covers needle tip automatically after use
- Looks and feels like a standard cannula

CODE	Dimension	Info	Pack Qty
121020	G24 0,7 x 19 mm	Yellow	50
121021	G22 0,9 x 25 mm	Blue	50
121022	G20 1,1 x 25 mm	Pink	50
121023	G20 1,1 x 33 mm	Pink	50
121024	G18 1,3 x 33 mm	Green	50
121025	G18 1,3 x 45 mm	Green	50
121026	G17 1,5 x 45 mm	White	50
121027	G16 1,7 x 50 mm	Gray	50
121028	G14 2,2 x 50 mm	Orange	50



Infusion needle Venflon Pro

- Improved insertion which increases the chance of first stick success, reducing the risk of thrombophlebitis
- The improved design of the injection port cap ensures easy opening and secure closing, whilst avoiding touch contamination of the injection port and valve
- Flexible wings stabilize the catheter in the vein and provide comfortable securement
- Catheters are compatible with BD Luer-Lok™ Syringes
- Sterilised by irradiation with a shelf life of 3 years
- CE marked and latex free
- Single use catheters



CODE	Dimension	Info	Pack Qty
130130	0,9 x 25 mm	Blue	50
130131	1,1 x 32 mm	Pink	50
130132	1,3 x 32 mm	Green	50
130133	1,3 x 45 mm	Green	50
130134	1,5 x 45 mm	White	50
130135	1,8 x 45 mm	Gray	50
130136	2,0 x 45 mm	Orange	50

Injection needles

- BD™ Microlance™ Stainless Steel Needles feature thin walls for optimal flow and fit Luer slip syringes
- Transparent with color code for fast and easier control of flow and size
- Sterile and individually wrapped for hygienic, single-use application
- Stainless steel needles with epoxy resin adhesives
- Do not contain: Latex, PVC or phthalates



CODE	Dimension	Info	Pack Qty
130200	0,3 x 13 mm	Yellow	100
130201	0,4 x 13 mm	Gray	100
130202	0,4 x 19 mm	Gray	100
130203	0,4 x 38 mm	Gray	100
121061	0,4 x 40 mm	Gray	100
130204	0,45 x 10 mm	Brown	100
130205	0,45 x 13 mm	Brown	100
130206	0,45 x 16 mm	Brown	100
130207	0,5 x 16 mm	Orange	100
130208	0,5 x 25 mm	Orange	100
130209	0,6 x 25 mm	Blue	100
130210	0,6 x 30 mm	Blue	100
130211	0,7 x 25 mm	Black	100
130212	0,7 x 30 mm	Black	100
130213	0,7 x 40 mm	Black	100
130214	0,7 x 50 mm	Black	100
130215	0,8 x 16 mm	Green	100
130216	0,8 x 25 mm	Green	100
130217	0,8 x 40 mm	Green	100
130218	0,8 x 50 mm	Green	100
121080	0,8 x 80 mm	Light green	100
130219	0,9 x 25 mm	Yellow	100
130220	0,9 x 40 mm	Yellow	100
130221	1,1 x 25 mm	Beige	100
130222	1,1 x 40 mm	Beige	100
130223	1,1 x 50 mm	Beige	100
130224	1,2 x 40 mm	Pink	100
130225	1,2 x 50 mm	Pink	100
130226	1,6 x 40 mm	White	100
130227	2,0 x 80 mm	Light green	100

Injection needle B.Braun Sterican

- Single use needle
- Minimal pain upon puncture
- Transparent Luer-Lock plastic hub
- Thin-walled needles
- Made of stainless, chromium-nickel steel
- Smooth surface with light silicone coating

CODE	Dimension	Info	Pack Qty
121081	0,8 x 120 mm	Green	100
121083	1,2 x 50 mm	Pink	100



Injection needle TSK LDS hub

- Aesthetic dermal filler cannulas with blunt tips have been developed for use with Hyaluronic Acid and non-Hyaluronic Acid dermal fillers
- The length and flexibility allow cosmetic practitioners to cover a larger surface area with a single-entry point
- Blunt tip dermal filler cannulas slide through subcutaneous layer rather than cutting through tissue and vessels like traditional aesthetic hypodermic needles

CODE	Dimension	Pack Qty
839906	33G 0,20 x 9mm	100
839907	33G 0,20 x 13mm	100



Injection needle TSK Steriglide

- Aesthetic dermal filler cannulas with blunt tips have been developed for use with Hyaluronic Acid and non-Hyaluronic Acid dermal fillers
- Steriglide creates less resistance than a traditional cannula tip. This reduces the cannula introduction force needed and the risk of bruising
- Highest quality steel with a proprietary surface treatment

CODE	Dimension	Pack Qty
839930	30G 0,30 x 25 mm	20
839931	27G 0,40 x 50 mm	20
839932	27G 0,40 x 38 mm	20
839933	27G 0,40 x 25 mm	20
839934	25G 0,50 x 50 mm	20
839935	25G 0,50 x 38 mm	20
839936	22G 0,70 x 70 mm	20
839937	22G 0,70 x 50 mm	20



Aspiration needle BD Nokor Admix

- BD filter needles feature a molded-in 5-micron filter membrane that reduces larger particles during fluid aspiration, reconstitution or transfer, and during medication withdraw from glass reservoirs

CODE	Dimension	Info	Pack Qty
130520	18G 1,25 x 40mm	Rosa	100
130521	16G 1,60 x 25mm	Lila	100

