


# Sterilization



### 3M Comply Indicator tape

- Indicator for hydrogen peroxide after sterilization cycle.
- The indicator change color after the product is exposed to sterilization.

CODE	Size	Pack Qty
100790	19 mm x 55 m	120



### Autoclave tape

CODE	Size	Info	Pack Qty
765020	19 mm x 50 m	With indicator	1
765021	25 mm x 50 m	With indicator	1
765010	19 mm x 50 m	Without indicator	1
765011	24 mm x 50 m	Without indicator	1



### Melacontrol Indicator (steam, time and temperature)


- Chemical indicator used as both a steam penetration test and to check the critical process parameters (steam, time and temperature) at each run in the autoclave.
- The package includes 250 indicators and new test body/box.

CODE	Pack Qty
107024	1 test body and 250 indicators



### Indicator (steam, time and temperature)

CODE	Pack Qty
107022	250 indicators



## Instrument basket

- Sterilization basket for instruments.
- Lid not included.

CODE	Size	Pack Qty
107029	210 x 150 x 40 mm	1
107031	310 x 210 x 40 mm	1
107032	420 x 310 x 40 mm	1



## Lid for instrument basket

- Lid for sterilization basket.

CODE	Size	Pack Qty
107039	210 x 150 mm	1
107041	310 x 210 mm	1
107042	420 x 310 mm	1



## Instrument cover

- Instrument cover in medical paper for steam or EtO sterilization.
- 15mm slits

CODE	Info	Pack Qty
510080	15 mm slits	1000



## Plastic instrument cover

- Transparent instrument cover in polypropylene.
- Autoclavable and adapted for plasma, ethylene oxide and steam sterilization.
- Disposable.

CODE	Info	Pack Qty
510085	15 mm slits	200
510086	30 mm slits	200



### Sterilization container

- Sterilization container intended for steam sterilization.
- Facilitates handling and storage of sterilized goods.
- Available in two sizes and can be stacked.
- Quick and easy replacement of paper filters in the lid and bottom.
- Made of aluminum.

CODE	Size	Pack Qty
107062	310 x 190 x 130 mm	1
107063	310 x 190 x 65 mm	1



### Moisture paper 80 g/m<sup>2</sup>

- Produced from pure cellulose.
- Can be safely used with sensitive instruments for steam sterilization.
- The high absorbency provides fast absorption of the condensed product, which facilitates the sterilization process.

CODE	Info	Pack Qty
510042	305 x 200 mm	200
510043	305 x 410 mm	200
510044	457 x 550 mm	200



### Sterilization wrap, Crepe paper green 60 g

CODE	Info	Pack Qty
150100	50 x 50 cm	500
150101	60 x 60 cm	500
150102	75 x 75 cm	250
150103	90 x 90 cm	250
150104	100 x 100 cm	250
150105	120 x 120 cm	100



## Halyard Quick Check Sterilization wrap

- Maintains sterility and aseptic presentation of the sterilized goods.
- Suggested usage based on wrap grade is given below.

CODE	Grade	Info	Suggested usage	Pack Qty
510230	H100	60 x 60 cm	Light linen packs/gowns, towel packs, basins	240
510231	H100	76 x 76 cm	Light linen packs/gowns, towel packs, basins	144
510232	H100	91 x 91 cm	Light linen packs/gowns, towel packs, basins	144
510233	H100	101 x 101 cm	Light linen packs/gowns, towel packs, basins	120
510234	H100	121 x 121 cm	Light linen packs/gowns, towel packs, basins	120
741040	H200	76 x 76 cm	Linen towels, gowns, light basin sets, single instruments	144
741041	H200	91 x 91 cm	Linen towels, gowns, light basin sets, single instruments	144
741042	H200	101 x 101 cm	Linen towels, gowns, light basin sets, single instruments	120
741043	H200	121 x 121 cm	Linen towels, gowns, light basin sets, single instruments	48
741044	H200	137 x 137 cm	Linen towels, gowns, light basin sets, single instruments	48



## Sterilization bag Cleantex NonWoven

- Manufactured in NonWoven material
- Used for heavy goods or instrument grids
- Disposable product, can be sterilized and opened once
- Used in the same way as a regular paper / plastic sterilization bag
- Easy to open
- Equipped with steam and EtO indicators



CODE	Info	Pack Qty
510320	210 x 420 mm	1000
510321	270 x 350 mm	1000
510322	270 x 450 mm	500
510323	320 x 500 mm	500
510324	420 x 600 mm	500

## Sterilization bag

- Transparent, multilayer film for easy identification.
- Reinforced film to avoid tearing when opening.
- Superior protection with 60g medical paper.
- Can be easily welded together with all bag welding machines.
- Strong sealing strength with visible weld seams.
- Marked with opening direction.
- Compliance with EN 868-5, ISO 11607.
- Disposable.



CODE	Info	Pack Qty
FP7515	75 x 150 mm	100
FP7520	75 x 210 mm	100
FP7527	75 x 280 mm	100
FP1020	100 x 200 mm	100
FP1025	100 x 240 mm	100
FP1027	100 x 280 mm	100
FP1030	100 x 300 mm	100
FP1040	100 x 420 mm	100
FP1520	150 x 200 mm	100
FP1527	150 x 280 mm	100
FP1530	150 x 320 mm	100
FP1535	150 x 360 mm	100
FP1540	150 x 400 mm	100
FP1552	150 x 520 mm	100
FP1660	160 x 600 mm	100
FP2027	200 x 320 mm	100
FP2040	200 x 420 mm	100
FP2550	250 x 500 mm	100
FP3050	300 x 500 mm	100
FP3057	300 x 580 mm	100

## Sterilization bag self-adhesive

- Self-adhesive bags do not require bag welding.
- Excellent adhesive, resistant to steam overexposure.
- Printed water-based non-toxic steam and EtO indicator.
- Folded sealing lip.

CODE	Info	Pack Qty
KP0720	75 x 200 mm	200
KP0923	90 x 230 mm	200
KP1024	100 x 240 mm	200
KP1325	135 x 255 mm	200
KP1428	140 x 280 mm	200
KP2042	200 x 330 mm	200
KP2540	250 x 400 mm	200

