

Suturing

Skin adhesive

Skin staplers

Sutures



Front Distribution

Medical Equipment

info@frontdistribution.com

www.frontdistribution.com

DermaBond wound glue Mini 0,36 ml

- DermaBond is a medical skin adhesive that is used to glue the sides of an incision or injury closed.
- It has shown to improve healing rates, reduce the chances of an infection, and the patient spends less time seeking medical attention.
- It can be used instead of stitches for small incisions and lacerations, but also in conjunction with deep dermal stitches.



CODE	Volume	Pack Qty
466000	0,36 ml	12

Skin adhesive Histoacryl

- Histoacryl® is available in two colours: translucent, special for facial application, and blue, which enables an easy control over the quantity applied in different indications approved.
- Closure of skin wounds without tension (including clean surgical incisions and incisions from minimally invasive surgery), and simple, thoroughly cleansed, trauma-induced lacerations.
- Sclerotherapy of large fundal varices.
- Fixation of hernia meshes, especially in inguinal hernia surgery.
- Wound closure in one minute and significant less pain than suture materials.
- Storage at room temperature



CODE	Volume	Info	Pack Qty
121248	0,5 ml	Transparent	5
121249	0,5 ml	Blue	5

Skin stapler Precise

- Quick and easy to use
- Elimination of hatching
- Easy removal for excellent cosmetic results
- Reduced incidence of infection
- Be able to select the number of staples required for each cut
- Reliable performance
- Choose between 5, 15 or 25 shaped staples per handle.

CODE	Info	Pack Qty
100660	DS5	12
100662	DS15	12
100661	DS25	12



Skin stapler Proximate Plus MD 55W

- Single use, sterile skin stapler
- 35 staples per device
- Staples have dimensions of 0,53 mm (diameter) x 3,9 mm (length) x 5,7mm (width)

CODE	Info	Pack Qty
PMR35	DS5	6



Staples remover BBraun

- Single use

CODE	Pack Qty
121235	6



Suture removal set

- One set contains:
 - 1 tweezers 127 mm
 - 3 dryers 30x30x30 mm
 - 1 suture knife
 - 1 tray

CODE	Pack Qty
701460	15



Ethibond Excel

- Braided Coated: For easy handling and secure knot tying. For smooth tie down and easy passage through tissue.
- Synthetic and non-absorbable: Provides permanent tensile strength retention in tissue.
- Indications: Ethibond sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurosurgical procedures. Typical areas of use include: heart valve fixation, ligament and tendon repair, retraction sutures.



CODE	Info	Pack Qty
EH7346H	Needle: 19 mm FS-2 Gauge: 3-0 Color: Green Length: 45 cm	36
EH7579H	Needle: 22 mm PB-1 Gauge: 3-0 Color: Green Length: 45 cm	36
JA6442H	Needle: 13 mm P-3 Gauge: 4-0 Color: White Length: 45 cm	36
EH7760H	Needle: 16 mm PS-3 Gauge: 4-0 Color: Green Length: 45 cm	36

Ethilon

- Monofilament: For smooth passage through tissue.
- Synthetic and non-absorbable: Provides prolonged tensile strength retention in tissue.
- Indications: Ethilon (Black and Clear) sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, microsurgical and neurosurgical procedures. Blue Ethilon sutures are intended for use in skin closure.
- Typical areas of use include: cuticular sutures, nerve adaptation, ophthalmology.

CODE	Info	Pack Qty
664H	Needle: 26 mm FS Gauge: 2-0 Color: Black Length: 45 cm	36
EH7826H	Needle: 24 mm FS-1 Gauge: 2-0 Color: Blue Length: 45 cm	36
1674H	Needle: 36 mm FSLX Gauge: 2-0 Color: Blue Length: 75 cm	36
628H	Needle: 60 mm KS Gauge: 2-0 Color: Black Length: 75 cm	36
3460TG	Needle: 76 mm LR Gauge: 2-0 Color: Black Length: 3x50 cm	12
663H	Needle: 24 mm FS-1 Gauge: 3-0 Color: Black Length: 45 cm	36
653H	Needle: 19 mm FS-2 Gauge: 3-0 Color: Black Length: 45 cm	36
1629H	Needle: 24 mm FS-1 Gauge: 4-0 Color: Black Length: 45 cm	-36
EH7144H	Needle: 19 mm FS-2 Gauge: 4-0 Color: Blue Length: 45 cm	36
662SLH	Needle: 19 mm FS-2S Gauge: 4-0 Color: Black Length: 45 cm	36
1162H	Needle: 16 mm FS-3 Gauge: 4-0 Color: Black Length: 45 cm	36
1894G	Needle: 19 mm PC-5 Gauge: 4-0 Color: Black Length: 45 cm	12



Monocryl Plus

CODE	Type	Info	Pack Qty
MCP4423H	Gauge 3-0	FS-2 - 70 cm White	36
MCP293H	Gauge 3-0	FS-2 - 45 cm White	36
MCP442H	Gauge 3-0	FS-1 - 70 cm White	36
MCP292H	Gauge 4-0	FS-2 - 45 cm White	36
MCP4422H	Gauge 4-0	FS-2 - 70 cm White	36
MCP494H	Gauge 4-0	P-3 - 70 cm White	36
MCP291H	Gauge 4-0	FS-1 - 70 cm White	36



Monocryl

CODE	Type	Info	Pack Qty
Y443H	Gauge 2-0	FS-1 - 70 cm White	36
Y213H	Gauge 5-0	RB-1 - 70 cm Purple	36
Y493H	Gauge 5-0	P-3 - 45 cm White	36
Y492H	Gauge 6-0	P-3 - 45 cm White	36



Stratafix Suture PGA-PCL

- Knotless tissue device that provides strong, secure closure appropriate for high-tension areas such as fascia.


CODE	Type	Info	Pack Qty
SXMD2B410	Gauge 3-0	FS-1 - 30 cm White	12



Stratafix Suture PDO

- Knotless tissue device that provides strong, secure closure appropriate for high-tension areas such as fascia.
- Absorberas mellan 120-180 dagar.

CODE	Type	Info	Pack Qty
SXPD2B411	0-0	MH - 30 cm Purple	12



Vicryl Plus

CODE	Type	Info	Pack Qty
VCP9361H	Gauge 0-0	V-30 - 70 cm Purple	36
VCP339H	Gauge 2-0	CT-1 - 70 cm	36
VCP442H	Gauge 3-0	FS-1 - 70 cm White	36



Vicryl Rapide

CODE	Type	Info	Pack Qty
V2930H	Gauge 3-0	FS-2 - 45 cm White	36
FH1642H	Gauge 4-0	FS-2 - 45 cm White	36
V2180H		SH-1 - 70 cm White	36
V2920H		FS-2S - 45 cm White	36
V2910H	Gauge 5-0	FS-2 - 45 cm White	36
V4930H		P-3 - 45 cm White	36
W9913	Gauge 6-0	P-1 - 45 cm White	12
V32H		P-1 - 45 cm Purple	36



Vicryl

CODE	Type	Info	Pack Qty
V443H	Gauge 2-0	FS-1 - 70 cm White	36
V453H		FS-1 - 70 cm Purple	36
V333H		CT-2 - 70 cm Purple	36
V602H		UR-6 - 70 cm Purple	36
V461H		X-1 - 70 cm Purple	36
V276H	Gauge 3-0	V-5 - 70 cm Purple	36
V311H		SH-1 - 70 cm Purple	36
V293H		FS-2 - 45 cm White	36
V393H		FS-2 - 45 cm Purple	36
V497H		PS-2 - 45 cm White	36
V416H		SH - 70 cm White	36
V332H		CT-2 - 70 cm Purple	36
V442H		FS-1 - 70 cm White	36
V423H		FS-2 - 70 cm White	36
V4416H		SH - 70 cm Black	36
V303H	Gauge 5-0	RB-1 - 70 cm Purple	36
V391H		FS-2 - 45 cm Purple	36
V421H		FS-2 - 70 cm White	36
V388H		FS-3 - 45 cm Purple	36
V490H		P-1 - 45 cm Purple	36
V493H		P-3 - 45 cm White	36
W9436		PC-4C - 45 cm White	36
PZ2571H	Gauge 6-0	RB-2 - 70 cm Purple	36
V489H		P-1 - 45 cm White	36
V570G		2x S-14 - 45 cm Purple	36

